

# Perfectionism

Being the pioneer of **folding mobility equipment**, FreeRider has been innovating powerful and beautiful mobility equipment to the riders.

Not only adopt the advantages of Luggie scooters also integrate the patented Coaxial **DMS tech** to improve the traditional power chair main problem in order to deliver superior **traction control management.** 





## Innovation but NOT RECREATION

Any products need to have innovation, but without documentation your innovation is ZERO.

Today, we're proud of Luggie chair is not belong to recreational product category so we keep our liability to ensure every consumers rights are protected under by high standard FDA requirement in the USA.

**Luggie Chair is FDA Approval Medical Devices: K 160139** 



## Armrest

Armrest height and width are adjustable for comfortable your position.

# Battery

Battery can be adopted by all Luggie lithium-ion freely. Options can be 8.5ah, 10.5ah and 16.5ah. Warranty are all covered 12 months.



## Caster

Casters are designed by quick-released. Will help you to put into trunk easily.

R

F

Ţ

ļ

0

ONE

Features

A pair of front anti-tip designing for power chair to deliver more safety and well-equipped with seat belt.

### TWO

Two seat height position to be adjusted

### THREE

Without any complicated steps and tools can be fold and expand in 3 second ready to go.



Stand-able foot platform & Quick-Release-able casters



### THINK AHEAD FOR YOU





P \_

=

=

)

T

!

N I

# UN 38.3

An International testing criteria for lithium battery or equipment containing lithium battery should undergo all testing criteria.

Our 8.5ah and 10.5ah lithium-ion battery and Luggie Chair itself are certified the UN 38.3 criteria. Battery weight is 4.5lbs.

Customers can be carry the light weight battery to airline.


It is **TSA approved.** 







Powered by 18650 Li-ion battery pack with certified IEC 62133 and UN 38.3. In addition, it adopts intelligent BMS battery management system which comes with various battery protections for the battery cell and your using safety.

# LUGGIE Li-ion Guide









#### **Medical Devices**



TravelRider is certified all the power wheelchair requirement, so it comes with K number based on FDA regulation. We submitted it as Medical Devices to ensure you using is safety and provide product liability insurance.

#### Efficacy



Adopted aircraft aluminium class as chassis and with zero-welding technique to make the strongest chassis than other power chair suppliers. It is smallest footprint and biggest weight capacity folding power chairs in the world.

#### **DMS System**



The patented and coaxial DMS system solve the traditional main problems. It's capable to delivers superior traction control management than others.

### Upgradable




We have different lithium-ion **battery options** to the meet your different requirements: 8.5ah, 10.5ah, and 16.5ah. You can choose bigger capacity if you need to have long journey.


#### Warranty



Lithium-ion battery is the most common as a portable power supply. However, we are the only one to provide **12 months warranty** instead of 6 months warranty from others.

Top Speed	3.6mph
Weight Cap	320lbs
Turning Radius	30''
Range	up to 22 mile (16.5ah)
Fully Charging	180 min
Overall Length	37''
Overall Width	23''
Caster	7'' solid PU
Rear Wheel	8" solid PU
Ground Clearance	2.36 (at motor)
Seat	18*16''
Total Weight	70lbs
Seat height adjustable	33" & 30"







2 0 1 9

# JOIN US

### FREERIDER USA

8696 Utica Ave Rancho Cucamonga, CA 91730

Email: sales@freeriderusa.com

TEL: 909-466-7996 FAX: 909-466-7990