

Q50R/carbon

Travel lighter



RIDICULOUSLY LIGHT: 32 LBS

One of the world's lightest carbon fiber folding power wheelchairs.

MODERN DESIGN

Sleek, ergonomic design with premium materials – the epitome of comfort and style. Available in 5 accent colors: Blue, Orange, Purple, Red, and Yellow.

GO UP TO 15 MILES

Bring an (optional) second battery along for the ride to double your range to up to 15 miles!

EFFORTLESS FOLDING

Flip the clip & push the backrest down. You don't even need to remove the battery.

COMPACT FOLDED SIZE

Folds to a compact size of 30.2" length x 22.4" width x 14.6" height, to fit small car trunks, buses, and closets.

SMOOTH RIDE

Its carbon fiber frame can naturally flex over uneven surfaces, providing unique shock-absorbing characteristics to deliver one of the smoothest, most efficient rides around.

LIGHT BUT INCREDIBLY STRONG

Can support a user weight of up to 300 lbs.

NARROW 22.4" WIDTH

Its slim width easily fits through narrow doorways and navigates busy rooms.

RELIABLE MOBILITY

250W high-efficiency brushless motors provide the power and reliability you need.



TECHNICAL DATA

DIMENSIONS	36.5" Long x 22.4" Wide
SEAT DIMENSIONS	17" Wide X 16.5" Deep
BACK HEIGHT	18.3"
MAX. RANGE	7.5 miles (one battery), 15 miles (with optional additional battery)
MAX. SPEED	3.7 MPH
BATTERY SIZE	10Ah lithium-ion
DRIVE WHEEL POSITION	Rear
WHEEL SIZE	8" Drive Wheel; 6.5" Casters
USER WEIGHT CAPACITY	300 lbs.
PRODUCT WEIGHT	Chair (without batteries): from 32 lbs., Battery (single): 3 lbs.

Refer to order form as specifications may vary based on configuration.

[Learn more](#)



SUNRISE
MEDICAL



QUICKIE